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DAKO 5140

This instruction manual is an integral part of the machine supply and must accompany it whenever it is moved. The manual must be kept carefully for the duration of the machine's life and must in all cases be available for

at least 10 years, it must therefore be kept in a known location and be made available to all the personnel concerned.
Do not connect or start up the machine before you have read this manual.


Fig.1


Fig.2

- 1 Main switch
- 2 Lid handle
- 3 Breakaway torque control knob
- 4 Portlight
- 5 Optical reader
- 6 Optical reader cable
- 7 Lid
- 8 Temperature setting device
- 9 Control panel
- 10 Side panel
- 11 Front panel
- 12 Dataplate
- 13 Power cable
- 14 Panel fastening screws
- 15 Water drain
- 16 Lock

- 1 Fixed counterweight
- 2 Mobile counterweight
- 3 Handwheel
- 4 Flask holder check pins
- 5 Flask holder
- 6 Crucible holding plate
- 7 Crucible
- 8 Flask
- 9 Crucible lid fastener
- 10 Flask bearing plate

TECHNICAL SPECIFICATIONS

The **DAKO 5140** machine is available in the following versions:

220 V 50 Hz	
220 V 60 Hz	
240 V 50 Hz	
Max. power rating	2860 W
Width	64cm + 5cm main cable
Depth	55 cm
Height	88cm + 5.5cm for optical reader
Net weight	89 Kg
Sound pressure level	<70 dB(A)

Time required for melting 30 gr of CrCo (four ingots weighing 10 g each) at the maximum power rating, with C15N crucible: approx. 55 sec.

PURPOSE OF THE EQUIPMENT

DAKO 5140 is an electronic casting machine equipped with microprocessor supervision system, and intended for casting the alloys to be used in the dental field. Its function is to induction melt the alloy in a suitable crucible and then, when the alloy is liquid, to inject it into a parallel mould, previously conditioned.

Any use other than above mentioned has to be considered as an improper use, which relieves the manufacturer from any responsibility for damages of any kind whatever, furthermore involving the user to forfeit with any guarantee right of the device.

PRECAUTIONS

- The machine must be used by properly trained adult persons.
- In consideration of the machine weight, any displacement must be carried out by means of a suitable mean.
- Avoid installing the machine close to a steam generator.
- During any displacement or transport, keep the machine upright (not horizontal). The machine must always stand on its feet.
- Wear safety glasses with dark lens (the glasses are included in the kit provided with the machine). The light produced by the melting metal could cause damages to the sight.
- Always use tongs to take off the crucible and the flask.
- Never start the casting process unless the crucible has been put in place. In case of use of a ceramic crucible, at least 15 gr of alloy must be herein.
- Avoid switching off the machine with heating coil up and hot crucible dipped in the heating coil.
- Check the turns of the heating coil are not in contact with each other.
- Make a flask large enough. The use of too small a flask could cause the early cooling of a part of the metal.
- Never use metal dirty of oil, grease or investment.
- Align the outlet of the crucible with the inlet of the flask, then slightly turn the crucible in the counter-clockwise direction of rotation.
- Make sure the electric system is equipped with ground safety cutout and protective earthing.
- Use only ceramic rods to touch the metal into the crucible during the casting process.

UNPACKING

The machine is delivered in a cardboard box laid on a wooden board. To remove the packaging, perform the following steps:

- Cut off the straps.
- Cut the adhesive tape and open the upper part of the cardboard box.
- Remove the staples from the base of the box.
- Remove the box slipping it upwards.
- Remove the board and place the machine on the floor.
- Do not try to open the lid. The lid can be opened only after the machine has been installed and plugged in. It is necessary to fill its water tank like explained at "charging with cooling water" and Plug the device into a suitable socket, like explained at "POWER SUPPLY"

POSITIONING

- The machine must be placed in a dry place free of dust, without any boilers or furnaces for flasks, in particular the gas ones especially if their draught is

not quite in order.

- Always place the machine on a flat surface. Check the surface flatness by means of a spirit level.

POWER SUPPLY

- ***The machine must be earthed by law.***
- ***Must be installed by a qualified electrician capable of checking that it is earthed.***
- ***The manufacturer declines all responsibility for any damage or injury caused by failure to comply with this requirement.***
- Plug the device into a suitable socket, having voltage and power rating corresponding to the rated output of the device, and fitted with ground.

ACTIVATION

Charging with cooling water

Take off the plug from the socket, if it's in. Remove the front panel (11, fig. 1), then pour approx. seven litres of distilled water into the tank. Replace the front panel, and then put the plug in the socket again.

Assembling the portlight

In case the portlight is not in place or was removed, follow this procedure:

1. Pull out from the box of the accessories the envelope containing the glass, the glass fastener, a grid and three fastening screws.
2. Place the glass over the porthole in contact with the lid, then place the grid and the glass fastener.
3. Secure the whole with the three fastening screws.

Assembling the centrifugation arm

If the arm has not been installed or it must be replaced, perform the following steps:

1. Pull the centrifugation arm out from the box of the accessories.
2. Switch on the machine by setting the main switch (1, fig. 5) to I.
3. Open the machine lid (7, fig. 1)
4. Undo completely the handwheel.
5. Insert the centrifugal arm and place on the driving shaft so that the lower slots rest on the drawbolt and the crucible side is turned towards the direction of the arrow printed on the driving shaft.
6. Put the two washers on the driving shaft, first the plain washer and then the indented one.
7. Tighten the central handwheel.

OPERATING INSTRUCTIONS

Please find below a detailed explanation of how to use the machine correctly. Afterwards, the standard working procedure is described, with a list of the steps that have to be carried out in succession.


Fig.3


Fig.4


Fig.5

1. ON/OFF switch.
2. Centrifuge Breakaway Torque Control
3. Melting Start key.
4. Coil Raise key.
5. Automatic Casting-Centrifugation key and centrifugation time setting.
6. Centrifugation Start key.
7. Lid locked ideogram.
8. Melting/Centrifugation Stop and Coil Lower key.
9. Programmable Value Increment key.
10. Programmable Value Decrement key.
11. FUNC key; no action takes place if this key is individually depressed.
12. REV key; no action takes place if this key is individually depressed.

Note: the simultaneous pressure of the FUNC and REV keys for 5 seconds gives the operator the possibility to access some configuration parameters and to change them. However, such change produces an alteration in the machine operation. Therefore, this combination of keys must never be used! Should this happen nevertheless due to an error, avoid pressing any key as, after a few seconds, the instrument is automatically reset to the normal operating mode.

13. ▼ key to decrement the temperature setting.
14. ▲ key to increment the temperature setting.
- A. Multifunction Led indicators:
 - melting power rating,
 - melting stabilisation time,
 - centrifugation time.
- B. Centrifugation time scale in seconds.
- C. Melting stabilisation time scale in seconds.
- D. Arm positioning Led indicator; it starts blinking when the centrifuge arm is correctly positioned, thus enabling the coil to raise. This led keeps lit with fixed light when the coil is up.
- E. Melting Led indicator: it keeps lit with fixed light during the melting process. It starts blinking when the lid opening breaks the melting process or a temperature control alarm is raised and when the melting time has to be entered.
- F. Automatic mode Led indicator: it keeps lit with fixed light during the entering of the related parameters or during the machine running in the Automatic mode. It blinks for five seconds before the centrifuge is started in the Automatic mode.
- G. Centrifugation Led indicator: it keeps lit with fixed light during the centrifuge operation. It starts blinking when the centrifugation time has to be entered.
- H. Lid locked Led indicator: it keeps lit with fixed light as long as the lid is held locked for safety reasons, i.e. when during the centrifuge arm rotation. It turns off after the centrifuge arm has come to a complete stop.
- I. STOP Led indicator: it's lit with fixed light when the machine is idle, i.e. when neither melting nor centrifugation is in progress.
- J. PV display (Left.): indicates the metal temperature in the crucible. A few seconds after the machine is switched on, the 800 °C temperature reading is displayed, that is. the lowest working temperature limit.
- K. SP display (Right.): indicates the temperature setting. It turns on when the machine is switched on and it displays the actual temperature setting.

The following indicator lights are installed on the instrument:

- AL: Temperature alarm indicator light. This lamp turns on when, because of a fault, the temperature to be stabilised is exceeded by approx. 50°C. In this case, heating is limited by the alarm itself.
- OUT1: Temperature control indicator light. It turns on when the machine is switched on. It blinks during the melting process in case the temperature to be stabilised is reached.
- °C: This lamp lit means that the temperature readings on the display are expressed in centigrade degrees.
- °F: This lamp lit means that the temperature readings on the display are expressed in Fahrenheit degrees.

Balancing of the centrifugation arm

- Lay a support (5, fig. 2) suitable to the flask to be used, on the centrifugation arm in order to get the support is firmly held between its locking pins (4, fig. 2), then put the flask therein (8, fig. 2).
 - Put the crucible with the metal into its seat (6, fig. 2).
 - Loosen the handwheel (3, fig. 2) of the arm.
 - Move the counterweight till the arm is balanced. The counterweight runs freely when it is pushed towards the central shaft. To move it in the opposite direction press the lever placed in its lower part.
 - Tighten the handwheel.
- Note: the arm where the mobile counterweight slides is marked with numbered positions to allow the balancing of all of flasks, by simply taking note of the counterweight location before the flask is entered to the heating furnace. On casting, the counterweight shall be moved to the position marked for the concerned flask.

Switching on

- Set the main switch (1, fig. 5) to I.
- The green light of the main switch turns on; all the led indicators of the keyboard light up for a few seconds; afterwards, the first six red Leds (A, fig. 4) and the STOP indicator light (H, fig. 4) keep lit and, provided the arm has already reached its position, also the indicator led corresponding to the Coil Raise push-button starts blinking (D, fig. 4).

Switching off

- To switch off the machine, press the main switch (1, fig. 7) to set it off. All the red Leds installed on the control panel and the green light of the main switch turn off.

If the heater coil is up and a hot crucible is dipped in the heating coil, first press the Stop push-button (8, fig. 4), and then switch off the machine.

Entering of the temperature value

To select the wanted temperature value, press either key 13 or 14 (see Fig. 5) according to whether you want to decrease or increase the temperature setting. To start a change in the temperature value, hold in the selected key for approx. two seconds. The temperature reading changes more quickly if the key is held in.

Note that the temperature entries must always be higher by some tenths of degrees than the maximum alloy casting range to get the metal is still fully melt when it reaches the flask inside.

Movements of the coil raise control lever

- From the bottom to the top (casting): rotate the arm and move it to the required position (the crucible must be over the heater coil) until the light of the coil raise push-button starts blinking (4, fig. 4), then depress the above mentioned push-button.
- From the top to the bottom: press the Stop push-button (8, fig. 4).

Melting start

To start the casting process, put the crucible on its support, then raise the coil and close the lid.

- To cast in the Manual mode, press the Casting key (3, fig. 4). Its Led indicator turns on.
- To cast in the Automatic mode, press the Auto key (5, fig. 4): the pilot light of the Casting Start key blinks. Afterwards, set the casting stabilisation time by means of the Increment or Decrement keys (9 and 10 on fig. 4) making reference to the scale C shown on fig. 4. Press the melting Start key (3, fig. 4): the pilot light of the Centrifugation key starts blinking. As above described, enter also the centrifugation time by making reference to the scale B shown on fig. 4. Finally, press the Centrifugation (6, fig. 4) and the Casting Start (3, fig. 4) keys.

At any time the operator may return to the Manual mode by pressing the Auto key (5, fig. 4).

The casting power rating may be changed, even if the casting process is already in progress, by means of the Increment or Decrement keys (9 and 10 on fig. 4).

The casting operation may start only on closed the heater coil is raised and the lid is closed. If the lid is opened during the casting operation, the process is automatically halted, then restarted when the lid is closed.

When the lid is opened the function AUTO is excluded and when the lid will be closed again the process will go on in manual mode

Never start the casting operation before the crucible is put on place. In case of use of a ceramic crucible, at least 15-g alloy must be therein.

Casting Stop

The casting process may be stopped in either way:

- press the Casting-Centrifugation Stop push-button (8, fig. 4);

- open the lid (when the lid is closed, the casting process is started again).

Centrifugation time

To set the centrifugation time, at the end of which the centrifuge is automatically stopped, press first the Auto key (5, fig. 4), and then the melting Start one (3, fig. 4): the light of the Centrifugation Start key starts blinking and the light of the Auto key keeps lit. Set the wanted time by operating the Increment and/or Decrement key (9 e 10, fig. 4), by making reference to the time scale B on figure 4. Finally, press the Centrifugation Start key (6, fig. 4) and the Auto one (5, fig. 4); their associated pilot light turns off.

Centrifuge Start

The centrifuge can be started only provided that all the safety conditions are met and the operating cycle is run correctly. Therefore, to start the centrifuge, follow this procedure:

1. in the Manual mode:

- with machine set in the melting mode and after closing the lid,
- press the Centrifugation key (6, fig. 4);

2. in the Auto mode (led of the Auto key (5, fig. 4) lit)

- the centrifuge starts immediately after the set temperature holding time has elapsed. Five seconds before the centrifugation, the Led of the Auto key (F, fig. 4) starts blinking and a sound signal is raised. To return to the Manual mode, simply press the Auto key.

During the centrifugation, the lid is locked and it cannot be opened (see also under "Lid opening").

Centrifuge Stop

The centrifuge may be stopped in either way:

- automatically: after a set time, the power supply to the centrifuge motor is cut off and the centrifugation arm coasts and stops;
- manually: to stop the arm before an automatic stop, press the Casting-Centrifugation Stop push-button (8, fig. 4).

Opening of the lid

The lid may be opened only on condition the Lid locked pilot light (H, fig. 4) is off and the machine is live. Under this condition, simply grip the handle and pull it up.

SELECTION OF THE CRUCIBLE

To get the best results from your casting machine, carefully apply the following guidelines.

Ceramic crucibles: C15N, C15R.

C15N: Standard crucible, standard mixture, recommended for casting of non-precious and CrCo alloys.

C15R: Ceramic crucible with liner suitable to insure longer life to the crucible. Especially recommended for non-precious alloys for ceramic. However, its performances are excellent also for other non-precious alloys.

The two above mentioned crucibles have to be pre-heated for at least 5 to 10 minutes before using them, as a rule in the same furnace where also the as-cast flasks are heated; don't forget to remove the slag left in the melting crucible. No vitrification is required.

Graphite crucibles (with graphite lining).

CG5: Standard graphite crucibles appropriate for casting of all the precious alloys suitable to be melt in such material.

CGV5: Sintered plumbago crucibles appropriate for casting of Palladium alloys and of all the precious alloys. Here all the precious-alloy castings result to be highly improved. When first using this crucible, care for heating it slowly; to this purpose, set the casting power rating at the lowest value (press the key 10 - fig. 4)) until - during the metal heating process - only one of the multifunction Led indicators is lit (A, fig. 4), taking care of keeping the temperature differential between heated metal and crucible to a minimum.

Note: the introduction of parts that, because of their shape and/or size, get stuck between the walls of the crucible may cause its breaking.

During the heating, these parts expand and exert a strong pressure against the walls of the crucible, thus causing its stress decay.

We remind you that you will automatically forfeit the right to the guarantee liding our equipment if you use thereon crucibles which are not of our production.

In order to avoid any dispute, all of our crucibles have been imprinted with our trademark.

WORKING PROCEDURES

In the Manual mode (Auto-key led off)

- Switch on the machine.
- Open the lid and balance the arm for the flask you intend to use.
- Rotate the arm until the light of the Coil Raise push-button starts blinking.
- Raise the coil by pressing the Coil Raise push-button.
- Adjust the starting acceleration of the centrifuge.
- Select the wished temperature.
- Close the lid.
- Set the melting power level.
- Press the melting key.
- Wait that the metal is melted.
- If this operation has not been carried out previously, open the lid and place the flask on the flask holder, in this case after the closing of the lid, it will be

necessary wait for the metal takes heat again.

- Press the centrifugation key.
- Let the centrifugal arm stops automatically once the programmed time is elapsed or press the stop key.
- Wait that the pilot lamp of locked lid goes off.
- Open the lid and take out the flask with suitable tongs.
- When all castings are finished switch off the machine by means of its main switch.

In the Automatic mode (Auto-key led on)

- Switch on the machine.
- Open the lid and balance the arm for the flask you intend to use.
- Rotate the arm until the light of the Coil Raise push-button starts blinking.
- Raise the coil by pressing the Coil Raise push-button.
- Adjust the starting acceleration of the centrifuge.
- Select the wished temperature.
- Place the flask on the flask holder
- Close the lid.
- Press the key AUTO its pilot lamp light up and the pilot lamp of the melting key blinks.
- Select the time of stabilization of the selected temperature.
- Press the melting key. The pilot lamp of the centrifugation key blinks.
- Select the centrifugal time.
- Press the centrifugation key, its pilot lamp goes off.
- set the melting power level.
- Press the melting key.
- Five seconds before centrifugation the pilot lamp of the AUTO key starts to blink.
- If the operator does not want to start automatically the centrifugation, press the key AUTO within this five seconds, in this way the machine return in not automatic functionment, otherwise the centrifuge will start automatically.
- Let the centrifugal arm stops automatically once the programmed time is elapsed or press the stop key.
- Wait that the pilot lamp of locked lid goes off.
- Open the lid and take out the flask with suitable tongs.
- When all castings are finished switch off the machine by means of its main switch.

sequence of sound signals. To identify the number of the Led that has turned on, start counting from the front of the machine backwards.

1. Led 1 lit: keyboard error. Call for assistance.
2. Led 2 lit: cooling water missing. Fill up with water, as previously explained. In case the water was just poured into, try to tilt the machine backwards and put a shim, 10 to 15 cm thick, under the front feet; leave the switched-on machine under this condition for approx. 10 minutes, then switch off the machine and put it down on the floor. Finally, switch on the machine again.
3. Led 1 and 2 lit: error of the coil lowering routine. Switch off the machine, then switch it on. Call the assistance service, if required.
4. Led 3 lit: lock open and lid closed. Switch off the machine, then switch it on again. Call the assistance service, if required.
5. Led 2 and 3 lit: lock closed and lid open. Switch off the machine, then switch it on again. Remember to close the lid before starting the centrifuge. Call for assistance, if required.
6. Led 4 lit: the maximum time set for the coil lifting motor has been exceeded. Start the machine again. Call for assistance, if required.
7. Led 1 and 4 lit: error of the coil raising routine. Switch off the machine, then switch it on again. Call the assistance service, if required.

- *The casting process is interrupted and a set of sound signals is raised.*

If several castings have been produced in succession, the cooling water temperature could be too high. Switch off the machine, then switch it on again after approx. one hour. If the trouble persists, call the assistance service.

- *The lid cannot be opened.*
The machine is not live. Power it up.

- *Vibrations of the machine during the centrifugation.*
The arm has not been correctly balanced. Balance it accurately according to the instructions given in this Manual.

FAULTS

- *The machine cannot be switched on.*
 1. Make sure the plug is in and the socket is live.
 2. Check the safety circuit breaker switch is set at I (ON). To this purpose:
 - remove the safety panel on the right (10, fig. 1).
 - set the magnetothermal switch to I.
 - replace the safety panel.
- The faults are identified by the microprocessor of the casting machine and displayed through the multifunction Led indicators (A, fig. 4) according to the following error codes, in conjunction with a

MAINTENANCE

- At least once in a year check the level of the distilled water. The level to be reached is marked on the water tank.
- Clean the glass of the portlight always before every casting process. If the glass is dirty, the optical reader cannot correctly take the temperature value.
- Once in a year remove the front panel in order to clean and lubricate with silicone lubricants the sliding surface of the coil raising system. The dust deposits could hinder the coil raising movement. Carefully clean the insulating plate that retains the two coil holding rods; finally put the front panel in place.

WARRANTY

Certificate of guarantee (approved by UNIDI)

1. Through the present document the manufacturer certifies a correct construction of the product and assures that first materials have been used; all necessary approvals and trials have been made so that the product is conform to the regulations in force. The guarantee on this product is valid for a period of 12 months starting on the day of delivery to the purchaser. The purchaser approves the guarantee regulations by filling in and signing his part and returning it to the manufacturer.

The guarantee only covers those parts to be substituted or replaced because of construction faults and doesn't cover labour expenses, travelling indemnity of technicians, transport expenses, packing etc.

The guarantee excludes damages or failures caused through a bad maintenance, a wrong feeding supply, negligences, unskillfulness or other causes not imputable to the manufacturer.

The guarantee excludes as well failures caused through omitted maintenance due to the negligences of the user.

The guarantee doesn't include any kind of refund versus persons or objects for damages, direct or indirect, due to eventual inefficiency of the equipment.

2. The guarantee forfeits automatically in case that the equipment has been repaired, modified or tampered by the purchaser himself or by unauthorized third parties.
3. For interventions under guarantee, the purchaser must address himself to the selling-agent, or the by him indicated technical-service or direct to the manufacturer. The guarantee gives the right to substitution of the defective part, but excludes the substitution of the entire equipment.
4. In case of contestation concerning the application of the guarantee regulations about the quality or conditions of the delivered equipments, the purchaser will in no way suspend or delay the payment of the price or the rates of price.
5. No refund can be requested by the purchases for the temporary stop of the equipment.
6. The guarantee forfeits if:
 - a) the equipment presents damages of dropping; exposure to flames; overturning of liquid; natural disasters; or any other causes not imputable to fabrication defaults;
 - b) the installation wasn't correct;
 - c) the connection to the main supply was wrong rated feeding voltage;
 - d) the serial no has been taken of, cancelled or altered.
7. The components, to be substituted under guarantee must be sent back to the manufacturer for

reimbursement following carefully his shipment instructions. For no reason we will accept the returning of equipments or spare parts of them without having given written authorization before. In case of non observation of these conditions the manufacturer reserves himself the right to send the goods back to the sender.


8. **Due to tax reasons the spare parts will be granted under guarantee only if we have received back the guarantee coupon properly filled in with all customer's data.**


Guarantee regulations


1. The non restitution of the guarantee certificate involves the immediate withdrawal of the certificate itself.

The guarantee doesn't cover labour expenses and travelling indemnity which will always be on charge of the purchaser.

2. The payment of the labour expenses, travelling indemnity, and right of appeal bill will be have to fulfilled immediate at the presentation of the bill itself. Not fulfilling the payment of the bill involves immediate withdrawal of the guarantee itself.
3. Neither the manufacturer or the dental deposit are obliged to provide substitutive equipment during the repair period.
4. For any other cases not foreseen in the present guarantee certificate and regulations reference is made to the italian civil code.

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Costruito a: **PORDENONE (PN) - ITALIA**
Made in:

costruito secondo le norme tecniche sotto riportate:
manufactured in conformity with the following standards:

Norme tecniche di riferimento: **EN 292-1, EN 292-2, EN 294, EN 60204-1, EN 50081-1, EN50082-1, EN 55011**
Technical norms of reference

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Direttive: **89/392CEE - 91/368CEE - 93/44CEE-93/68CEE-89/336CEE-92/31CEE-73/23CEE**
Directives:

20/05/02

Addi - Date

Bruno Bortolus

Amministratore Delegato – Managing Director

NOTE

This image shows a full page of white paper with horizontal dotted lines. The lines are evenly spaced and run across the width of the page, providing a guide for handwriting practice. There are no margins, text, or other markings on the page.

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SARATOGA S.P.A.

Manufacturing of dental surgeries and laboratories equipments

via Malignani, 14 - 33107 Pordenone - ITALY

Tel. +39.0434.572600 r.a. Fax +39.0434.572477

[http: www.saratoga-on-line.com](http://www.saratoga-on-line.com) e-mail: info@saratoga-on-line.com